

Lipogro

Cholesterol – Lipoprotein Complex, 0.1 micron filtered; 1L in PETG media bottle

Product number: RMB-001-001

Lot Number: 20120321G

Date of Manufacture: March 2012

Expiration date: March 2014

Storage temperature: 2 – 8°C

Not for human use

Basis for Release

Test / Method	Specification	Result
Cholesterol Concentration	9.00 – 11.00 g/L	9.84 g/L
Total Protein	13.0-18.0 g/L	14 g/L
pH	7.0 – 8.4	8.1
Sodium	≤ 50 mEq/L	Below Range
Potassium	0 – 10 mEq/L	Below Range
Chloride	≤ 6mEq/dL	Below Range
Calcium	≤ 5.0 mg/dL	<5 mg/dL
Endotoxin	< 50 EU/mL	1.4 EU/mL
Sterility	No growth	Pass
Virus Testing per 9CFR 113.53	None Detected	Pass
Virus Testing-Porcine Parvovirus	None Detected	Pass
Appearance	Clear amber solution, visibly free of particulates	Pass

Approval

Quality Assurance – Print Name: <i>Emily Villegas</i>	Quality Assurance – Sign Name: <i>Emily Villegas</i>	Date: <i>24 Apr 12</i>
--	---	---------------------------

Lipogro is derived from bovine serum, processed from blood collected from healthy animals, sourced from GBR 1-rated countries.

Characterization assays are performed on various lots of this product, and will be provided upon request.