

Bovine Serum Albumin, lyophilized powder

Product number: RMB-002-101 Lot Number: 20120410-S

Date of Manufacture: 04/2012 Expiration date: 04/2017

Storage temperature: 2 – 8°C

For evaluation purposes only – not for human or animal consumption

Basis for Release

Test / Method <i>Tests performed in 4% solution unless otherwise noted</i>	Result
pH	7.4
Solubility	pass
Protein, Total	3.8 g/dL
Purity	98%
Calcium	<6 mg/dL
Chloride	<25 mEq/L
Potassium	<1.0 mEq/L
Sodium	< 20 mEq/L
Ash	1.5%
Microbial	0 cfu/mL
Moisture	1.0%

Approval		
Quality Assurance – Print Name: <i>Abbigail Miller</i>	Quality Assurance – Sign Name: <i>Abbigail J. Miller</i>	Date: <i>4 Oct 12</i>

Bovine albumin is derived from bovine serum; processed from blood collected from healthy animals, sourced from GBR-1 rated countries. Characterization assays are performed on various lots of this product, and will be provided upon request.