

Equine Serum Albumin

Part Number: RMB-008-001 Lot Number: 20130613EA
 Date of Manufacture: 06/2013 Date of Expiration: 06/2018
 Storage temperature: 2 – 8°C

For research or manufacturing only. Not for human or animal consumption

Basis for Release

Test / Method <i>Tests performed in 4% solution unless otherwise noted</i>	Result
Origin	US Source
Appearance	Off White Crystal Powder
pH (1% solution)	7.36
Total Protein	3.3 g/dL
Calcium	< 5.0 mg/dL
Potassium	< 1.0 mmol/L
Sodium	13 mmol/L
Chloride	< 10 mmol/L
Purity	94%

Approval		
Quality Assurance – Print Name: <i>Emily Villegas</i>	Quality Assurance – Sign Name: <i>Emily Villegas</i>	Date: <i>06/26/13</i>

Equine albumin is derived from equine serum, processed from blood collected at a USDA licensed establishment and donor animals were found to be healthy at the time of collection.