

Maxichol®, Bovine Cholesterol Concentrate; 0.1µm filtered

Product number: MXL-BDG-XXX Lot Number: 20130726D
 Date of Manufacture: 07/2013 Expiration date: 07/2015
 Storage temperature: 2 – 8°C

For further manufacturing only. Not for human or animal consumption.

Basis for Release

Test / Method	Result
Cholesterol Concentration	1529 mg/dL
Total Protein	3.0 g/dL
pH	8.1
Sodium	22 mmol/L
Potassium	< 1.0 mmol/L
Chloride	<10 mmol/L
Calcium	< 5.0 mg/dL
Appearance	Clear amber solution, free of visible particulates

Approval

Quality Assurance – Print Name: <i>Emily Villeges</i>	Quality Assurance – Sign Name: <i>[Signature]</i>	Date: <i>09/12/13</i>
--	--	--------------------------

Maxichol® is derived from bovine serum, processed from blood collected from healthy animals, sourced from the US.