

### Bovine Gamma Globulin

Product number: BGG-BBZ-XXX      Lot Number: 20150424IG  
 Date of Manufacture: 04/2015      Expiration date: 04/2020  
 Storage temperature: 2 – 8°C



*Not for human or animal consumption*

#### Basis for Release

Test / Method <i>Tests performed in 5% solution unless otherwise noted</i>	Result
Origin	United States
Appearance	Off White
pH	6.79
Purity	98%
Calcium	< 2.0 mg/dL
Chloride	< 25 mEq/L
Potassium	< 0.2 mEq/L
Sodium	< 20 mEq/L
Microbial	Passed
Moisture	2.2%
Total Protein	4.5 g/dL

#### Approval

Quality Assurance – Print Name: <i>Tern Sorenson</i>	Quality Assurance – Sign Name: <i>TC Sorenson</i>	Date: <i>08/25/15</i>
---	--	--------------------------

Bovine Gamma Globulin is derived from bovine serum, processed from blood collected at a USDA regulated abattoir. The animals received ante and post mortem inspection at the abattoir and were found free of contagious diseases.