

Bovine Serum Albumin, Fatty Acid Poor

Product number: RMB-002-203

Lot Number: 51019ALD

Date of Manufacture: 05/2010

Expiration date: 05/2015

Storage temperature: 2 – 8°C

Not for human or animal use

Basis for Release

Test / Method <i>Tests performed in 4% solution unless otherwise noted</i>	Result
Appearance	Light powder
pH (1% solution)	6.75
Purity	98.1%
Calcium	< 5.0 mg/dL
Chloride	< 25.0 mEq/L
Potassium	< 0.2 mEq/L
Sodium	26.8 mEq/L
Fatty Acid	0.45 mg/g
Microbial	Passed
Moisture	1.1%
Filtration	Solution sterile filtered through 0.2µm prior to lyophilization

Approval		
Quality Assurance – Print Name: <i>Emily Villegas</i>	Quality Assurance – Sign Name: <i>Emily Villegas</i>	Date: <i>14 May 12</i>

Bovine albumin, Fatty Acid Poor, is derived from bovine serum, processed from blood collected from healthy animals. Characterization assays are performed on various lots of this product, and will be provided upon request.